

3514700 Clino® CR1 SC-GMP Cleaning trolley for Cleanroom

**Individualization
on request!**



Product description

Cleaning trolley with dosing and impregnation system MopScoop GMP and disposal for highly sensitive areas made of stainless steel.

- Cleaning method: dosing and impregnation system MopScoop GMP
- No service life problems thanks to just-in-time impregnation
- Touch-free operating concept
- Made of cleanroom-suitable stainless steel, all components autoclavable
- With disposal unit



ISO	3	4	5	6	7	8	9
GMP			A/B		C	D	

Product specifications

Article name Article number	Clino® CR1 SC-GMP 3514700
Working dimension	78 x 53 x 70 cm
Resting dimension	78 x 53 x 70 cm
Basic components	Basic module (stainless steel suitable for cleanroom use) High-end finish Autoclavable Special rollers 100 mm
Workspace	1 MopScoop dosing unit including stainless steel lid 1 MopTen stainless steel container with MopJump system 1 device holder
Disposal area	Unit 1 x 35 L
Sterility	The sterility of the trolley depends on the sterility of its individual parts. Our sales department will be happy to help you!
Care instructions	The cleaning trolley should be completely cleaned after use and stored in a dry place
Packing unit Selling unit	1 Pieces 1 PU
Quantity in carton	1 Pieces
Carton size Carton weight	78 x 61 x 106 cm 15 kg
Pallet	80 x 123 x 182 cm
Customs tariff number	87168000
Country of origin	Germany
EAN-Number	4065237109069
Compatible accessories	Our sales department will be happy to advise you on suitable accessories

Industry Sectors

Assembly Hygiene



Diese Klassifizierung ist eine unverbindliche Empfehlung. Der Anwender muss die individuelle Eignung und Verträglichkeit für seinen individuellen Anwendungsbereich prüfen. Der Zustand des Materials ist regelmäßig zu überprüfen, beschädigte Produkte sind umgehend auszutauschen. Alle Angaben ohne Gewähr. Änderungen und Abweichungen vorbehalten.

3514700 Clino® CR1 SC-GMP

Cleaning trolley for Cleanroom

1x	3530171	Module CR 1/2 - Facelift
1x	3530142	Bag support for CR4
4x	3500100	Stainless steel swivel castor 100 mm
1x	3530187	Stem for MopTen and system box
1x	3500022	MopTen
1x	3500210	MopJump for system box EasyMop GMP and MopTen
1x	2017031	Device holder 22 mm
1x	3500740	Trigger for MopScoop GMP
1x	3500710	MopScoop GMP container
1x	3500782	Lid for MopScoop GMP
1x	3500756	Seal for MopScoop GMP dosing tube
1x	3500722	Holding frame for MopScoop GMP
2x	3500735	Storage tray for MopScoop GMP tray
1x	3500753	MopScoop GMP dosing flask 200 ml
2x	3500755	Seal for MopScoop GMP dosing piston
1x	3500750	MopScoop GMP dosing tube
1x	3500762	Nut for MopScoop GMP
1x	3500730	Tray for MopScoop GMP
1x	0899120	Work instruction
1x	0899100	Cardboard small

Diese Klassifizierung ist eine unverbindliche Empfehlung. Der Anwender muss die individuelle Eignung und Verträglichkeit für seinen individuellen Anwendungsbereich prüfen. Der Zustand des Materials ist regelmäßig zu überprüfen, beschädigte Produkte sind umgehend auszutauschen. Alle Angaben ohne Gewähr. Änderungen und Abweichungen vorbehalten.

PFENNIG Reinigungstechnik GmbH | Heubachstraße 1 | 87471 Durach | Germany
www.pfennig-reinigungstechnik.com | info@pfennig.reinigungstechnik.com | +49 (0) 831 561 22 0

Status: 25.10.2023

